

Analysis Report

Official laboratory analysis summary for the submitted sample and associated quality-control review.

SAMPLE IGF-1 LR3	RECEIVED DATE May 18, 2026	ANALYSIS DATE May 22, 2026	REPORT GENERATED May 25, 2026
STRENGTH 1mg	MANUFACTURER PepticoresAminos	BATCH NUMBER PC-IGF1-0526E	LAB CODE 673-1
CLIENT www.pepticoresaminos.us			

SAMPLE INFORMATION

IGF-1 LR3**1MG**FORM **White lyophilized powder in a glass vial**SAMPLE SUBMISSION **Sample provided by customer**BATCH **PC-IGF1-0526E**CAP / CRIMP COLOR **Transparent/blue**RECEIVED DATE **May 18, 2026**

SAMPLE IMAGE


ANALYTICAL SUMMARY

IDENTITY	IGF-1 LR3
PURITY	99.441%
QUANTITY	1.04mg
BATCH	PC-IGF1-0526E
MANUFACTURER	PepticoresAminos

RP-HPLC-UV CHROMATOGRAM (214 NM)

Detection: UV 214 nm | Runtime: 20.0 min

Sample ID: IGF-1 LR3
Report ID: 2024-02-01-0001
Method: IGF-1 LR3 RP-HPLC-UV
Detector: UV 214 nm | Runtime: 20.0 min



METHOD

TIME	WATER + 0.1% TFA	ACETONITRILE + 0.1% TFA
0min	95%	5%
2min	95%	5%
10min	40%	60%
12min	5%	95%
20min	5%	95%

TECHNICAL NOTE

Analysis performed by RP-HPLC-UV at 214 nm using a TFA water/acetonitrile gradient adapted from IGF-1 LR3 peptide screening conditions. Total runtime: 20.0 minutes.

COMMENTS

The sample remains within the defined analytical specifications for identity, purity and impurity profile based on the applied RP-HPLC-UV method.

ANALYSIS & METHODOLOGY

IGF-1 LR3 RP-HPLC-UV ANALYSIS

HPLC-UV: High-Performance Liquid Chromatography with Ultraviolet detection. Additional testing includes bioburden screening, endotoxin analysis and heavy metals assessment.

BIOBURDEN

TEST	RESULT	UNIT	REPORTING LIMIT
Total Aerobic Microbial Count USP <61>/Eur. Ph. 2.6.12. Plate Count Method	Not detected	CFU/g	>= 1000
Total Yeast and Mold Count USP <61>/Eur. Ph. 2.6.12. Plate Count Method	Not detected	CFU/g	>= 100

ENDOTOXIN ANALYSIS

TEST	RESULT	UNIT	REPORTING LIMIT
Bacterial Endotoxin USP<85>/ Eur. Ph. 2.6.14. Bacterial Endotoxin Chromogenic Test	< 0.001	EU/mg	> 0.5

HEAVY METALS

TEST	RESULT	UNIT	REPORTING LIMIT
Arsenic Elemental Impurities Screening	Not detected	ppm	>= 1.5
Cadmium Elemental Impurities Screening	Not detected	ppm	>= 0.5
Cobalt Elemental Impurities Screening	Not detected	ppm	>= 25
Lead Elemental Impurities Screening	Not detected	ppm	>= 1.5
Nickel Elemental Impurities Screening	Not detected	ppm	>= 25
Quicksilver Elemental Impurities Screening	Not detected	ppm	>= 1.5
Vanadium Elemental Impurities Screening	Not detected	ppm	>= 25

TECHNICAL APPENDIX

This appendix documents the analytical methodology, instrumentation and acceptance criteria applied for the evaluation of the sample.

METHOD SPECIFICATION

PARAMETER	IGF-1 LR3 RP-HPLC-UV METHOD
ANALYTICAL MODE	RP-HPLC-UV peptide identity and purity screen for IGF-1 LR3
COLUMN	Phenomenex Biozen Peptide Polar C18, 150 x 2.1 mm, 3 um equivalent
MOBILE PHASE A	Water + 0.1% TFA
MOBILE PHASE B	Acetonitrile + 0.1% TFA
FLOW RATE	0.4 mL/min
DETECTION	UV 214 nm
INJECTION VOLUME	0.5 uL
COLUMN TEMPERATURE	60 C
RUNTIME	20.0 min

PARAMETER	IGF-1 LR3 RP-HPLC-UV METHOD
SAMPLE DILUENT	Aqueous organic diluent compatible with peptide RP-HPLC screening
SAMPLE PREPARATION	Diluted, mixed and clarified before chromatographic review

INSTRUMENT PLATFORM

PARAMETER	IGF-1 LR3 RP-HPLC-UV PLATFORM
SYSTEM TYPE	High-resolution LC platform
DETECTOR	UV/VIS detector configured at 214 nm
ACQUISITION	Chromatographic acquisition and integration software
REVIEW MODE	Large peptide retention-time and peak purity profile review
WORKFLOW NOTE	Used for IGF-1 LR3 RP-HPLC-UV identity and purity screening

ANALYTICAL CRITERIA

PARAMETER	ACCEPTANCE FRAMEWORK	BASIS
IDENTITY	Retention-time and profile agreement with reference expectations	Chromatographic identity review
PURITY	NLT 98.0% unless a stricter report-specific specification is declared	Integrated RP-HPLC-UV purity profile
QUANTITY	Measured content reviewed against the declared sample strength	Report-level analytical summary
BIOBURDEN	Not detected or within stated reporting limits	Microbial screening table
ENDOTOXIN	Below stated reporting limit / internal screening threshold	Endotoxin analysis table
HEAVY METALS	Below individual reporting limits where screened	Elemental impurities screening table

This appendix documents the analytical methodology, instrumentation and acceptance criteria applied for the evaluation of the sample.

VERIFICATION

Verify this report through the official Synaptica Analytics verification page using the details below.

Verification URL synaptica-labs.com/verify-report

Report ID SYN-2026-004961

Verification Key VK-JMKY-2L7X

SCAN TO VERIFY



DIGITAL SIGNATURE



DIGITALLY SIGNED BY:

Martin Saar

Date: 2026.05.25

16:20:27 +02'00'

Director

info@synaptica-labs.com

Analysis date: May 22, 2026

Report generated: May 25, 2026

Analytical testing performed by Synaptica Analytics -

Analytical Services Division

Synaptica Analytics
SYN-2026-004961
Laboratory Analysis Report
VK-JMKY-2L7X

VERIFY AT
synaptica-labs.com/verify-report